

## Content of Modify the AutoChk Timeout Value.vbs ( Site 1 )

' Description: Sets the auto-delay time for Autochk.exe to 30 seconds.

```
strComputer = "."
Set objWMIService = GetObject("winmgmts:"
    & "(impersonationLevel=impersonate)!\\" & strComputer & "\root\cimv2")

Set colAutoChkSettings = objWMIService.ExecQuery _
    ("Select * from Win32_AutochkSetting")

For Each objAutoChkSetting in colAutoChkSettings
    objAutoChkSetting.UserInputDelay = 30
objAutoChkSetting.Put_
Next
```