

Content of WMI MethodLoadUnload Query.js (Site 1)

```
var wbemFlagReturnImmediately = 0x10;
var wbemFlagForwardOnly = 0x20;

var arrComputers = new Array("");
for (i = 0; i < arrComputers.length; i++) {
    WScript.Echo();
    WScript.Echo("=====");
    WScript.Echo("Computer: " + arrComputers[i]);
    WScript.Echo("=====");

    var objWMIService = GetObject("winmgmts:\\\\" + arrComputers[i] + "\\root\\WMI");
    var collItems = objWMIService.ExecQuery("SELECT * FROM MethodLoadUnload", "WQL",
        wbemFlagReturnImmediately | wbemFlagForwardOnly);

    var enumItems = new Enumerator(collItems);
    for (; !enumItems.atEnd(); enumItems.moveNext()) {
        var objItem = enumItems.item();

        try { WScript.Echo("EventGuid: " + (objItem.EventGuid.toArray()).join(", ")); }
        catch(e) { WScript.Echo("EventGuid: null"); }
        WScript.Echo("EventSize: " + objItem.EventSize);
        WScript.Echo("EventType: " + objItem.EventType);
        WScript.Echo("Instancelid: " + objItem.Instancelid);
        WScript.Echo("KernelTime: " + objItem.KernelTime);
        WScript.Echo("MethodFlags: " + objItem.MethodFlags);
        WScript.Echo("MethodIdentifier: " + objItem.MethodIdentifier);
        WScript.Echo("MethodSize: " + objItem.MethodSize);
        WScript.Echo("MethodStartAddress: " + objItem.MethodStartAddress);
        WScript.Echo("MethodToken: " + objItem.MethodToken);
        WScript.Echo("ModuleID: " + objItem.ModuleID);
        WScript.Echo("MofData: " + objItem.MofData);
        WScript.Echo("MofLength: " + objItem.MofLength);
        try { WScript.Echo("ParentGuid: " + (objItem.ParentGuid.toArray()).join(", ")); }
        catch(e) { WScript.Echo("ParentGuid: null"); }
        WScript.Echo("ParentInstancelid: " + objItem.ParentInstancelid);
        WScript.Echo("ReservedHeaderField: " + objItem.ReservedHeaderField);
        WScript.Echo("ThreadId: " + objItem.ThreadId);
        WScript.Echo("TimeStamp: " + objItem.TimeStamp);
        WScript.Echo("TraceLevel: " + objItem.TraceLevel);
        WScript.Echo("TraceVersion: " + objItem.TraceVersion);
        WScript.Echo("UserTime: " + objItem.UserTime);
    }
}
```