

Content of WMI Win32_PerfFormattedData_NETFramework_NETCLRJit Query.vbs (Site 1)

```
On Error Resume Next
strComputer = "."
Set objWMIService = GetObject("winmgmts:\\" & strComputer & "\root\cimv2")
Set collItems = objWMIService.ExecQuery("Select * from Win32_PerfFormattedData_NETFramework_NETCLRJit",,48)
For Each objItem in collItems
    Wscript.Echo "Caption: " & objItem.Caption
    Wscript.Echo "Description: " & objItem.Description
    Wscript.Echo "Frequency_Object: " & objItem.Frequency_Object
    Wscript.Echo "Frequency_PerfTime: " & objItem.Frequency_PerfTime
    Wscript.Echo "Frequency_Sys100NS: " & objItem.Frequency_Sys100NS
    Wscript.Echo "ILBytesJittedPersec: " & objItem.ILBytesJittedPersec
    Wscript.Echo "Name: " & objItem.Name
    Wscript.Echo "NumberofILBytesJitted: " & objItem.NumberofILBytesJitted
    Wscript.Echo "NumberofMethodsJitted: " & objItem.NumberofMethodsJitted
    Wscript.Echo "PercentTimeinJit: " & objItem.PercentTimeinJit
    Wscript.Echo "StandardJitFailures: " & objItem.StandardJitFailures
    Wscript.Echo "Timestamp_Object: " & objItem.Timestamp_Object
    Wscript.Echo "Timestamp_PerfTime: " & objItem.Timestamp_PerfTime
    Wscript.Echo "Timestamp_Sys100NS: " & objItem.Timestamp_Sys100NS
    Wscript.Echo "TotalNumberofILBytesJitted: " & objItem.TotalNumberofILBytesJitted
Next
```